

SCALABLE AI-DRIVEN IMAGE MANAGEMENT SOLUTION

EndoVault's versatile platform is designed to meet the imaging needs of a single room center to a large enterprise hospital system. EndoVault gives immediate access to images & videos, and is a vendor-neutral solution. EndoVault is a DICOM viewer application that integrates with all major EHR, VNA & PACS, eliminating manual scanning and dependence on paper or film.



SAVING MONEY WITH A MULTI-SPECIALTY IMAGING SOLUTION

Scalable software creates a cost-effective solution for an ASC, outpatient center, or large enterprise hospital system.

ENDOVAULT IMAGE MANAGEMENT APPLICATION

Can be utilized to meet the needs for multi-specialities including:

- ✓ Gastroenterology
- ✓ Pulmonology
- ✓ ENT
- ✓ Oncology
- ✓ Pain Management
- ✓ Urology
- ✓ General Surgery
- ✓ Cardiology
- ✓ Orthopedics
- ✓ Pathology
- ✓ Obstetrics/Gynecology
- ✓ Plastic Surgery
- ✓ Podiatry
- ✓ Dermatology
- ✓ Infusion Therapy

The EndoVault infrastructure shared across multiple specialties reduces the number of vendors to manage and support, while saving cost to implement and maintain.

ENHANCED EFFICIENCY & SEAMLESS SOLUTIONS

Intelligently capture and label images during a procedure using your voice and Argus AI. Capture from up to four sources simultaneously, view up to 1080p HD images in BMP, JPEG, TIF, and/or DICOM format, capture & view video in up to 4K in AVI format. From any scope, capture using HDMI, HD/SDI, DVI, S-video, RGB, AV and/or composite video sources. Compare and annotate images and video. DICOM visible light, worklist, and storage compliant that supports encapsulated PDF. IHE certified XDS.b-i & XDS.b cross enterprise document sharing.



+1 518-831-8000



sales@endosoft.com



www.endosoft.com

ENDOVAULT® INTEGRATED PACS

EndoVault® PACS provides accelerated image retrieval, economical storage and simultaneous access to images across multiple healthcare sites.

- ✓ Built-in DICOM viewer gives providers real-time access to images from PACS
- ✓ Supports DICOM Visible Light still image and streaming video capture.
- ✓ Supports DICOM encapsulated PDF



MODALITY COMPATIBILITY

Compatible with all modalities (CT, MRI, Nuclear Medicine, Ultrasound, Endoscopy, X-Ray, Angiography, Mammography, Computed Radiology, Microscopes, Digital & Analog Camera, etc.) all manufacturers and all DICOM 3.0 and NEMA 2 file types.

CERTIFICATIONS & COMPLIANCES

EndoVault® holds multiple certifications & compliances for all needs including ONC Health IT Certified Electronic Health Record (EHR) Application, IHE Compliant, Vista CPRS, Vista Imaging, NEDi2 Compliant, Infoway Certified Edition 2017, and National Cancer Screening Register (NCSR) Certified.



+1 518-831-8000



sales@endosoft.com



www.endosoft.com